

# Salmon Quantification - Comprehensive Reference Guide

## Overview

Salmon is a fast, alignment-free tool for quantifying transcript and gene expression from RNA-seq data. It uses a k-mer-based quasi-mapping approach combined with an Expectation-Maximization (EM) algorithm to estimate transcript abundances without requiring full genome alignment.

**Website:** <https://combine-lab.github.io/salmon/>

**Publication:** Patro et al. Nature Methods 2017

**Current Version:** 1.10.0+

**Applicable to:** RNA-seq, single-cell RNA-seq, metatranscriptomics, spatial transcriptomics

---

## What Salmon Does

### Core Functions

1. **Quasi-mapping** - Maps reads to transcripts using k-mer matching (not full alignment)
2. **Transcript Quantification** - Estimates abundance for each transcript (TPM and counts)
3. **Gene-level Aggregation** - Can aggregate to gene level (or use tximport)
4. **Multi-mapper Resolution** - Probabilistically assigns reads mapping to multiple transcripts
5. **Bias Correction** - Corrects for GC, sequence, and positional biases

### Salmon vs Traditional Methods

Feature	Salmon	STAR + featureCounts	kallisto
<b>Speed</b>	Very fast (2-5 min)	Slow (20-40 min)	Very fast (3-8 min)
<b>Accuracy</b>	Excellent	Excellent	Excellent
<b>Output</b>	Transcript counts	Gene counts	Transcript counts
<b>BAM file</b>	No	Yes	No
<b>Memory</b>	8-12GB	30-50GB	4-8GB
<b>Bias correction</b>	Built-in	Manual	Limited
<b>Multi-mappers</b>	EM algorithm	Discard or count once	EM algorithm
<b>Gene families</b>	Handles well	Problematic	Handles well

## Salmon Quantification Approach

### K-mer Based Quasi-Mapping

**Traditional alignment** (STAR, HISAT2):

Read → Find seeds → Extend alignment → Score → Report best alignment

Time: ~10-30 minutes per sample

Output: Full BAM file with exact positions

**Salmon quasi-mapping:**

Read → Hash k-mers → Find transcript matches → Score compatibility → Report

Time: ~2-5 minutes per sample

Output: Transcript abundance estimates (no BAM)

### How It Works

Step 1: K-mer Extraction

Read: ATGCGTACGATCGATCG...

K-mers (k=31): ATGCGTACGATCGATCGATCGATCGATC

TGCGTACGATCGATCGATCGATCGATCGATCG

GCGTACGATCGATCGATCGATCGATCGATCGA

...

Step 2: Transcript Matching

K-mer → Index lookup → Find transcripts containing this k-mer

Build compatibility graph: Read ↔ {Transcript1, Transcript2, ...}

Step 3: Probabilistic Assignment (EM Algorithm)

Initialize: Uniform probability across compatible transcripts

Iterate:

1. E-step: Assign reads to transcripts weighted by current estimates

2. M-step: Update transcript abundances based on assignments

Until convergence

Result: Probability distribution of each read across transcripts

→ Fractional counts for each transcript

### Why EM Algorithm?

**Problem:** Multi-mapping reads

Read maps to:

- Transcript A (gene X)
- Transcript B (gene X, different isoform)
- Transcript C (pseudogene of X)

Traditional methods:

Discard (lose information)  
Count once (underestimate)  
Count all (overestimate)

Salmon's EM:

Assign fractionally based on:  
- Transcript abundance  
- Sequence compatibility  
- Fragment length distribution  
Iterate until stable  
Maintains total read count

**Example:**

Initial:

Transcript A: 1000 reads (100 unique, 900 ambiguous)  
Transcript B: 100 reads (10 unique, 90 ambiguous)

EM realizes:

- A is 10x more abundant  
- Ambiguous reads more likely from A

After EM:

Transcript A: ~945 reads (probabilistic assignment)  
Transcript B: ~55 reads (probabilistic assignment)

---

## Salmon Quantification Command

### Basic Command

```
salmon quant \  
  --index salmon_index \  
  --libType A \  
  --mates1 R1.fastq.gz \  
  --mates2 R2.fastq.gz \  
  --output sample_quant \  
  --threads 4
```

### Command Breakdown

salmon quant \	<i># Quantification mode</i>
--index salmon_index \	<i># Path to Salmon index</i>
--libType A \	<i># Auto-detect library type</i>
--mates1 R1.fastq.gz \	<i># Forward reads (PE)</i>
--mates2 R2.fastq.gz \	<i># Reverse reads (PE)</i>

```
--output sample_quant \      # Output directory
--threads 4                  # Number of CPU cores
```

Single-end data:

```
salmon quant \
  --index salmon_index \
  --libType A \
  --unmatedReads reads.fastq.gz \
  --output sample_quant \
  --threads 4
```

## Key Parameters Explained

### Library Type (--libType)

**Purpose:** Specifies strand orientation of reads

**Auto-detection (RECOMMENDED):**

```
--libType A    # Auto-detect library type
```

Salmon will examine first 1-10 million reads and determine:

- **SE Libraries:** SF, SR, or U
- **PE Libraries:** ISF, ISR, or IU

**Library type codes:**

Code	Description	Read Orientation
<b>SF</b>	Single-end Forward	Read matches transcript strand
<b>SR</b>	Single-end Reverse	Read is reverse complement
<b>U</b>	Unstranded	Read could be either orientation
<b>ISF</b>	Inward Stranded Forward	R2 matches transcript, R1 is reverse
<b>ISR</b>	Inward Stranded Reverse	R1 matches transcript, R2 is reverse
<b>IU</b>	Inward Unstranded	Reads could be either orientation

**Common protocols:**

Protocol	Library Type
<b>Illumina TruSeq</b>	ISR (reverse stranded)
<b>dUTP method</b>	ISR (reverse stranded)
<b>Ligation-based</b>	ISR or ISF
<b>SMARTer</b>	Varies (use auto)
<b>Nextera</b>	Varies (use auto)
<b>Old protocols</b>	IU (unstranded)

### Verification:

```
# Check detected library type
cat sample_quant/lib_format_counts.json
```

```
# Example output
{
  "expected_format": "ISR",
  "compatible_fragment_ratio": 0.9823,
  "num_compatible_fragments": 9823000,
  "num_assigned_fragments": 10000000,
  ...
}
```

### If auto-detection fails:

```
# Manually specify (not recommended)
--libType ISR    # If you know it's reverse-stranded PE
```

---

### Selective Alignment (`--validateMappings`)

**Purpose:** More stringent mapping validation

**How it works:**

**Without `--validateMappings`** (default quasi-mapping):

1. Find k-mer matches
  2. Count matches per transcript
  3. Assign based on counts
- Fast but may include spurious matches

**With `--validateMappings`** (selective alignment):

1. Find k-mer matches
  2. Build alignment chains
  3. Score alignment quality
  4. Filter low-quality matches
  5. Assign based on validated matches
- ~20% slower but more accurate

**Accuracy improvement:** ~2-5% more accurate quantification

### Recommendation:

- Always use for publication-quality results
- Can skip for quick exploratory analysis

```
salmon quant \
  --validateMappings \    # RECOMMENDED for accuracy
```

```
--index salmon_index \  
...
```

---

## Bias Corrections

Salmon provides three types of bias correction:

**1. GC Bias (--gcBias)** **Problem:** PCR amplification favors certain GC content

Without correction:

High GC transcripts: Under-estimated  
Low GC transcripts: Over-estimated

Example:

GC-rich gene (70% GC): 500 reads → Actually 600 reads  
GC-poor gene (30% GC): 500 reads → Actually 400 reads

**How Salmon corrects:**

1. Model GC content vs fragment count
2. Learn bias curve from data
3. Re-weight fragments based on GC content
4. Output corrected abundances

**Command:**

```
--gcBias    # Enable GC bias correction
```

**Impact:** ~3-5% improvement in accuracy

**When to use:**

- PCR-amplified libraries
  - Non-UMI protocols
  - Publication-quality analysis
  - UMI-based methods (already corrected)
- 

**2. Sequence Bias (--seqBias)** **Problem:** Random hexamer priming isn't truly random

Random hexamer priming:

Primers: NNNNNN (supposedly random)  
Reality: Some hexamers bind better than others

Bias example:

GCGCGC: Binds strongly → Over-represented

ATATAT: Binds weakly → Under-represented

Result:

Transcripts starting with preferred hexamers → Over-estimated  
Transcripts starting with poor hexamers → Under-estimated

How Salmon corrects:

1. Learn hexamer preferences from read starts
2. Model position-specific bias
3. Correct abundances based on sequence context

Command:

```
--seqBias    # Enable sequence-specific bias correction
```

Impact: ~2-3% improvement in accuracy

When to use:

- Random hexamer priming protocols
  - Any RT-PCR based method
  - Oligo-dT priming (different bias pattern)
- 

**3. Positional Bias (--posBias)** **Problem:** Fragments not uniformly distributed along transcripts

Causes:

- 3' bias: RNA degradation, poly-A selection
- 5' bias: Some library prep methods
- Coverage holes: Secondary structure

Without correction:

3' biased sample:

Gene A (5' region expressed): Under-estimated  
Gene B (3' region expressed): Over-estimated

How Salmon corrects:

1. Model fragment position distribution
2. Learn coverage pattern along transcripts
3. Re-weight fragments based on position

Command:

```
--posBias    # Enable positional bias correction
```

Impact: ~1-2% improvement (varies by sample quality)

When to use:

- Degraded RNA samples
- Poly-A selected libraries (mild 3' bias)
- High-quality, uniform coverage

---

**Combined Bias Correction**    **Standard practice:** Enable all three

```
salmon quant \
  --validateMappings \
  --gcBias \
  --seqBias \
  --posBias \
  --index salmon_index \
  --libType A \
  --mates1 R1.fastq.gz \
  --mates2 R2.fastq.gz \
  --output sample_quant \
  --threads 4
```

**Total improvement:** ~5-10% more accurate quantification

**Trade-off:** ~10-15% slower runtime (still very fast overall)

---

### Mapping Parameters

**--minAssignedFragments**    **Purpose:** Minimum fragments assigned to a transcript for it to be reported

**Default:** 10

**Why it matters:**

Low-count transcripts (1-9 fragments):

- High uncertainty
- Likely noise or mapping artifacts
- Unstable estimates

Filtering helps:

- Reduces false positives
- Improves statistical power
- Cleaner downstream analysis

**Recommendations:**

```
# Standard analysis
--minAssignedFragments 10    # Default, good balance

# Sensitive analysis (detect low-abundance)
```



```
--minAssignedFrag 1      # Keep everything

# Conservative analysis
--minAssignedFrag 50     # Only well-supported transcripts
```

---

**--rangeFactorizationBins** Purpose: Controls EM algorithm precision

**Default:** 4

**What it does:**

Bins for fragment length distribution:

More bins → More precise → Slower

Fewer bins → Less precise → Faster

Default (4): Good balance

Higher (8-16): Marginal improvement, longer runtime

Lower (2): Faster but less accurate

**Recommendation:** Keep default unless you have specific needs

---

## Output Files

Salmon creates an output directory with multiple files:

### Directory Structure

```
sample_quant/
  quant.sf                # Main output: Transcript abundances
  quant.genes.sf          # Gene-level abundances (if --geneMap used)
  lib_format_counts.json  # Library type detection results
  meta_info.json          # Run metadata and parameters
  cmd_info.json           # Command line executed
  aux_info/
    ambig_info.tsv        # Ambiguously mapped fragment info
    expected_bias.gz      # Learned bias models
    fld.gz                # Fragment length distribution
    meta_info.json        # Auxiliary metadata
    observed_bias.gz      # Observed bias patterns
  logs/
    salmon_quant.log       # Detailed log file
```

---

### quant.sf (Main Output)

**Format:** Tab-separated values (TSV)

**Columns:**

Column	Name	Description
1	<b>Name</b>	Transcript ID (e.g., ENST00000000001)
2	<b>Length</b>	Transcript length in base pairs
3	<b>EffectiveLength</b>	Effective length accounting for biases
4	<b>TPM</b>	Transcripts Per Million (normalized abundance)
5	<b>NumReads</b>	Estimated number of reads from this transcript

**Example:**

Name	Length	EffectiveLength	TPM	NumReads
ENST00000000001.1	2000	1850	1250.50	10000.5
ENST00000000002.1	1500	1350	500.25	3500.2
ENST00000000003.1	3000	2850	2500.75	25000.8

**Column details:**

#### Length

- Raw transcript length from FASTA
- Includes all exons (spliced)
- Fixed value for each transcript

#### EffectiveLength

- Adjusted length accounting for:
  - Fragment length distribution
  - Positional bias
  - Edge effects (fragments can't start/end anywhere)

Formula:  $\text{EffectiveLength} = \text{Length} - \text{mean fragment length} + 1$   
where  $\text{mean fragment length}$

**Example:**

Transcript length: 2000bp  
Mean fragment length: 200bp  
Effective length:  $2000 - 200 + 1 = 1801\text{bp}$

Why: A 200bp fragment can't start in last 199bp of transcript

## TPM (Transcripts Per Million)

- Normalized abundance measure
- Comparable across samples
- Sums to 1 million within each sample

Formula:

1. Counts per base:  $\text{NumReads} / \text{EffectiveLength}$
2. Normalize:  $(\text{Counts per base} / \text{Sum of all counts per base}) \times 1,000,000$

Properties:

- Sum of all TPMs = 1,000,000
- Accounts for transcript length
- Comparable across samples
- NOT raw counts (don't use for DESeq2/edgeR directly)

## TPM interpretation:

TPM = 1000:	1 out of every 1000 transcripts in the sample
TPM = 100:	Moderate expression
TPM = 10:	Low expression
TPM = 1:	Very low expression
TPM < 1:	Barely detectable

## NumReads

- Estimated number of reads from transcript
- Can be fractional (due to multi-mapping)
- Sum across transcripts = total mapped reads

Example:

Read maps to Transcript A and B

If A is 3x more abundant than B:

Transcript A: +0.75 reads

Transcript B: +0.25 reads

Result: NumReads can be non-integer

---

## lib\_format\_counts.json

**Purpose:** Library type detection results

**Example:**

```
{  
  "expected_format": "ISR",  
  "compatible_fragment_ratio": 0.9823,  
}
```

```

    "num_compatible_fragments": 9823000,
    "num_assigned_fragments": 10000000,
    "num_consistent_mappings": 9823000,
    "num_inconsistent_mappings": 177000,
    "MSF": 50000,
    "MSR": 127000,
    "MU": 0,
    "SF": 25000,
    "SR": 25000,
    "U": 0,
    "ISF": 88500,
    "ISR": 9823000,
    "IU": 0
  }

```

#### Key fields:

- `expected_format`: Detected library type
- `compatible_fragment_ratio`: Fraction of reads matching expected orientation
- High ratio ( $>0.95$ ): Confident detection
- Low ratio ( $<0.80$ ): Mixed or poor quality

#### Red flags:

- Compatible ratio  $<0.80$ : Library type unclear
- Multiple types with similar counts: Mixed libraries
- Unexpected type: Wrong protocol or sequencing issue

---

#### `meta_info.json`

**Purpose:** Run metadata and parameters

#### **Example:**

```

{
  "salmon_version": "1.10.0",
  "samp_type": "paired",
  "opt_type": "vbem",
  "quant_errors": [],
  "num_libraries": 1,
  "library_types": ["ISR"],
  "frag_dist_length": 1000,
  "num_targets": 234567,
  "num_bootstraps": 0,
  "mapping_type": "mapping",
  "index_seq_hash": "a1b2c3d4e5f6...",

```

```

    "index_name_hash": "f6e5d4c3b2a1...",
    "index_decoy_seq_hash": "1a2b3c4d5e6f...",
    "index_keep_duplicates": false,
    "keep_duplicates": false
}

```

#### Important fields:

- **salmon\_version**: Version used (important for reproducibility)
- **num\_targets**: Number of transcripts quantified
- **library\_types**: Detected library type(s)
- **num\_bootstraps**: If bootstrapping was run

---

#### aux\_info/ Directory

Contains auxiliary information for downstream tools:

##### fld.gz

- Fragment length distribution
- Used by tximport for length correction
- Histogram of insert sizes

##### expected\_bias.gz / observed\_bias.gz

- Learned bias models
- Used for bias correction
- Can be visualized for QC

##### ambig\_info.tsv

- Information about ambiguously mapped fragments
- Useful for debugging low mapping rates

---

## Downstream Analysis with tximport

### Why Use tximport?

**Problem:** Salmon outputs transcript-level abundances, but we often want gene-level

#### Solutions:

1. Sum transcript counts directly → Loses statistical information
2. Use longest isoform only → Ignores real biology
3. Use tximport → Properly aggregates with statistical awareness

## tximport Advantages

- Aggregates transcripts to genes correctly
  - Preserves inferential uncertainty
  - Accounts for transcript length differences
  - Compatible with DESeq2/edgeR
  - Handles multi-isoform genes properly
- 

## Basic tximport Usage

### Step 1: Prepare tx2gene mapping

```
library(tximport)
library(dplyr)

# Option A: From Ensembl BioMart
library(biomaRt)
mart <- useEnsembl("ensembl", dataset = "hsapiens_gene_ensembl")
tx2gene <- getBM(
  attributes = c("ensembl_transcript_id", "ensembl_gene_id"),
  mart = mart
)
colnames(tx2gene) <- c("transcript_id", "gene_id")

# Option B: From GTF file
library(GenomicFeatures)
txdb <- makeTxDbFromGFF("annotation.gtf")
k <- keys(txdb, keytype = "TXNAME")
tx2gene <- select(txdb, k, "GENEID", "TXNAME")
colnames(tx2gene) <- c("transcript_id", "gene_id")

# Option C: Manual (if you have a file)
tx2gene <- read.table("tx2gene.txt", header = TRUE)
```

### Format:

transcript_id	gene_id
ENST000000000001	ENSG000000000001
ENST000000000002	ENSG000000000001
ENST000000000003	ENSG000000000002

---

### Step 2: Collect Salmon output files

```
# List all quant.sf files
samples <- c("sample1", "sample2", "sample3", ...)
```

```
files <- file.path("salmon_output", samples, "quant.sf")
names(files) <- samples

# Verify files exist
all(file.exists(files))
```

---

### Step 3: Import with tximport

```
# Import transcript abundances
txi <- tximport(
  files,
  type = "salmon",
  tx2gene = tx2gene,
  countsFromAbundance = "lengthScaledTPM" # Recommended
)

# Output structure
names(txi)
# [1] "abundance"           # Gene-level TPMs
# [2] "counts"              # Gene-level estimated counts
# [3] "length"              # Average transcript length per gene
# [4] "countsFromAbundance" # Method used
```

**countsFromAbundance options:**

Option	Description	Use Case
"no"	Use estimated counts directly	Standard
"scaledTPM"	Scale TPMs to library size	Simple normalization
"lengthScaledTPM"	Scale TPMs & adjust for length	<b>RECOMMENDED</b>
"dtuScaledTPM"	For DTU analysis	Isoform switching

---

### Step 4: Create DESeq2 object

```
library(DESeq2)

# Sample metadata
coldata <- data.frame(
  sample = samples,
  condition = c("control", "control", "treated", ...),
  row.names = samples
)
```

```

# Create DESeqDataSet from tximport
dds <- DESeqDataSetFromTximport(
  txi,
  colData = coldata,
  design = ~ condition
)

# Run differential expression
dds <- DESeq(dds)
results <- results(dds)

```

---

## Step 5: Explore results

```

# View results
head(results)

# MA plot
plotMA(results)

# PCA
vsd <- vst(dds)
plotPCA(vsd, intgroup = "condition")

# Top genes
res_ordered <- results[order(results$padj), ]
head(res_ordered, 20)

```

---

## Advanced tximport Features

### Import with bootstraps (if run)

```

# Run Salmon with bootstraps
salmon quant \
  --numBootstraps 30 \    # 30 bootstrap samples
  ...

# Import with uncertainty
txi <- tximport(
  files,
  type = "salmon",
  tx2gene = tx2gene,
  txOut = FALSE,          # Gene-level
  varReduce = TRUE        # Compute inferential variance
)

```



```
# Now txi includes variance estimates
names(txi)
# Includes "variance" element
```

**When to use bootstraps:**

- Publication-quality DEG analysis
  - When uncertainty quantification matters
  - Low-count genes
  - Exploratory analysis (overkill)
- 

## Transcript-level analysis

```
# Keep transcript-level (don't aggregate to genes)
txi_transcripts <- tximport(
  files,
  type = "salmon",
  txOut = TRUE,           # Keep transcript-level
  countsFromAbundance = "lengthScaledTPM"
)

# Use for isoform-level analysis
dds_tx <- DESeqDataSetFromTximport(
  txi_transcripts,
  colData = coldata,
  design = ~ condition
)
```

**Use cases:**

- Isoform switching analysis
  - Alternative splicing detection
  - Transcript-specific regulation
- 

## Interpreting Salmon Logs

### Mapping Rate

Check salmon\_quant.log:

```
[info] Mapping rate = 85.23%
```

**Interpretation:**

Mapping Rate	Quality	Action
>80%	Excellent	Proceed
70-80%	Good	OK, but check contamination
60-70%	Acceptable	Investigate low mapping
<60%	Poor	Check index, contamination, degradation

#### Common causes of low mapping:

- Wrong index (different species or annotation version)
- Contamination (bacterial, adapter, rRNA)
- Degraded RNA
- Poor sequencing quality

#### Fragment Length Distribution

Check aux\_info/fld.gz:

```
# Visualize in R
library(dplyr)
library(ggplot2)

fld <- read.table(gzfile("sample_quant/aux_info/fld.gz"))
colnames(fld) <- c("length", "count")

ggplot(fld, aes(x = length, y = count)) +
  geom_line() +
  labs(title = "Fragment Length Distribution",
       x = "Fragment Length (bp)",
       y = "Count")
```

#### Expected patterns:

Normal distribution (good):

100 200 300 (bp)

Bimodal (mixed libraries):

100 200 300

Very narrow (size selection issue):

200

Very wide (poor quality):

50      300      500

---

## Library Type Consistency

Check across samples:

```
# Extract library types from all samples
for sample in sample1 sample2 sample3; do
    echo -n "${sample}: "
    grep "expected_format" ${sample}/lib_format_counts.json
done
```

```
# Expected: All samples same type
sample1: "expected_format": "ISR"
sample2: "expected_format": "ISR"
sample3: "expected_format": "ISR"
```

```
# Red flag: Mixed types
sample1: "expected_format": "ISR"
sample2: "expected_format": "IU"      # ← Problem!
sample3: "expected_format": "ISR"
```

If inconsistent:

- Check library prep protocol
  - Verify not mixing different library types
  - May indicate sample swap or mislabeling
- 

## Resource Requirements

Memory (RAM)

---

Component	Requirement
<b>Salmon index in RAM</b>	5-8GB (human)
<b>Quantification buffer</b>	2-4GB

Component	Requirement
<b>Total</b>	8-12GB

#### Scaling by species:

- Mouse: 6-10GB
- Zebrafish: 4-6GB
- Drosophila: 2-4GB
- C. elegans: 1-2GB
- Arabidopsis: 2-3GB

#### CPU Cores

##### Scaling efficiency:

- 1 core: Baseline
- 4 cores: 3.5x faster (good efficiency)
- 8 cores: 6x faster (diminishing returns)
- 16 cores: 8x faster (poor efficiency)

**Recommendation:** 4-8 cores optimal

#### Disk Space

##### Per sample:

- Input FASTQ: 2-8GB (compressed)
- Output directory: 50-200MB
- Temporary files: Minimal

**For 100 samples:** ~5-20GB total output

#### Time

Sample Size	4 cores	8 cores
<b>10M reads</b>	2 min	1.5 min
<b>25M reads</b>	4 min	2.5 min
<b>50M reads</b>	7 min	4 min
<b>100M reads</b>	12 min	7 min

**With bias correction:** +10-20% time

## Workflow-Specific Considerations

### Bulk RNA-seq

#### Standard settings:

```
salmon quant \  
  --validateMappings \  
  --gcBias \  
  --seqBias \  
  --posBias \  
  --libType A \  
  --index salmon_index \  
  --mates1 R1.fastq.gz \  
  --mates2 R2.fastq.gz \  
  --threads 4 \  
  --output sample_quant
```

#### Expected metrics:

- Mapping rate: >75%
  - Library type: Consistent across samples
  - Fragment length: 150-300bp mean
- 

### Single-cell RNA-seq

#### Special considerations:

```
# alevin mode (for UMI-based scRNA-seq)  
salmon alevin \  
  --index salmon_index \  
  --libType ISR \  
  --mates1 R1.fastq.gz \  
  --mates2 R2.fastq.gz \  
  --chromiumV3 \  
  --output sample_alevin \  
  --threads 8
```

*# Cell barcode + UMI*  
*# cDNA read*  
*# Or other chemistry*

#### Differences:

- Uses alevin mode (not quant)
  - Handles UMIs and cell barcodes
  - Outputs count matrix (genes  $\times$  cells)
  - No bias correction (UMIs handle that)
-

## Metatranscriptomics

### Challenges:

- Multiple species in index
- High multi-mapping
- Taxonomic assignment needed

### Recommendations:

```
salmon quant \  
  --validateMappings \  
  --allowDovetail \  
  --recoverOrphans \  
  --libType A \  
  --index metatranscriptome_index \  
  --mates1 R1.fastq.gz \  
  --mates2 R2.fastq.gz \  
  --threads 8 \  
  --output sample_quant
```

### Downstream:

- Use taxonomic classification
  - Aggregate by species/genus
  - Account for shared genes
- 

## Troubleshooting

### Low Mapping Rate (<50%)

#### Diagnostic steps:

##### 1. Check index matches data:

```
# Verify species  
head -n 1 salmon_index/info.json
```

```
# Check number of transcripts  
grep "num_targets" sample_quant/meta_info.json
```

##### 2. Check for contamination:

```
# Run FastQC  
fastqc R1.fastq.gz R2.fastq.gz
```

```
# Check adapter content  
# Check overrepresented sequences
```

##### 3. Test with subset:

```

# Take first 100K reads
zcat R1.fastq.gz | head -400000 | gzip > R1.subset.fq.gz
zcat R2.fastq.gz | head -400000 | gzip > R2.subset.fq.gz

# Run Salmon
salmon quant --index salmon_index --libType A \
  --mates1 R1.subset.fq.gz --mates2 R2.subset.fq.gz \
  --output test_quant

# Check mapping rate
cat test_quant/logs/salmon_quant.log | grep "Mapping rate"

```

---

## Library Type Detection Fails

Symptoms:

```

{
  "expected_format": "IU",
  "compatible_fragment_ratio": 0.45,    // Low!
  ...
}

```

Causes & solutions:

1. **Truly unstranded library:**
    - OK if protocol is unstranded
    - Check library prep method
  2. **Mixed orientation:**
    - Sample contamination or mixing
    - Check sample sheet
  3. **Very low depth:**
    - Not enough reads for confident detection
    - Use `--libType ISR` (or appropriate type) manually
  4. **Poor quality:**
    - High error rate confuses detection
    - Check FastQC results
- 

## Out of Memory

**Error:** "std::bad\_alloc" or "Cannot allocate memory"

**Solutions:**

1. **Increase RAM allocation:**

- Allocate 16GB minimum
- 32GB for very large transcriptomes

2. **Reduce index size** (rebuild with --sparse):

```
salmon index \  
  --sparse \                # Sparse index (less memory)  
  --transcripts gentrome.fa \  
  --index salmon_index_sparse
```

3. Use --reduceGCMemory:

```
salmon quant \  
  --reduceGCMemory \        # Reduce garbage collection overhead  
  --index salmon_index \  
  ...
```

---

## Very Low Fragment Count

**Symptoms:** Most transcripts have <10 fragments assigned

**Causes:**

1. **Very low sequencing depth:**

- Check total read count
- May need more sequencing

2. **Wrong index:**

- Using transcriptome index when should use genome
- Or vice versa

3. **Severe degradation:**

- Fragments too short
  - Check fragment length distribution
- 

## Inconsistent Results Across Samples

**Check:**

1. **Library types match:**

```
grep "expected_format" */lib_format_counts.json
```

2. **Mapping rates similar:**

```
grep "Mapping rate" */logs/salmon_quant.log
```



3. Fragment lengths similar:

```
# Compare fld.gz across samples
```

4. Same index used:

```
grep "index_seq_hash" */meta_info.json | sort -u
```

---

## Best Practices

### Standard RNA-seq Quantification

```
salmon quant \  
  --validateMappings \  
  --gcBias \  
  --seqBias \  
  --posBias \  
  --libType A \  
  --index salmon_index \  
  --mates1 R1.fastq.gz \  
  --mates2 R2.fastq.gz \  
  --threads 4 \  
  --output sample_quant
```

---

### Quick QC Run (Speed Priority)

```
salmon quant \  
  --libType A \  
  --index salmon_index \  
  --mates1 R1.fastq.gz \  
  --mates2 R2.fastq.gz \  
  --threads 8 \  
  --output sample_quant
```

*# No bias corrections → 2x faster*

---

### Publication-Quality (Accuracy Priority)

```
salmon quant \  
  --validateMappings \  
  --gcBias \  
  --seqBias \  
  --posBias \  
  --numBootstraps 30 \  
  --libType A
```

*# Uncertainty quantification*

```
--index salmon_index \  
--mates1 R1.fastq.gz \  
--mates2 R2.fastq.gz \  
--threads 8 \  
--output sample_quant
```

---

## Related Documentation

- **Salmon Index:** docs/salmon\_index.md - Building the index
  - **tximport:** <https://bioconductor.org/packages/tximport/> - Importing to R
  - **DESeq2:** docs/deseq2.md - Differential expression
  - **Salmon Manual:** <https://salmon.readthedocs.io/>
- 

**Document Version:** 1.0

**Last Updated:** January 2026

**Salmon Version:** 1.10.0+

**Applicable to:** Bulk RNA-seq, scRNA-seq (alevin mode), metatranscriptomics