Package 'oligo'

May 22, 2024

Version 1.68.0

Title Preprocessing tools for oligonucleotide arrays

Author Benilton Carvalho and Rafael Irizarry

Contributors Ben Bolstad, Vincent Carey, Wolfgang Huber, Harris Jaffee, Jim MacDonald, Matt Settles, Guido Hooiveld

Maintainer Benilton Carvalho <benilton@unicamp.br>

Depends R (>= 3.2.0), BiocGenerics (>= 0.13.11), oligoClasses (>= 1.29.6), Biobase (>= 2.27.3), Biostrings (>= 2.35.12)

Imports affyio (>= 1.35.0), affxparser (>= 1.39.4), DBI (>= 0.3.1), ff, graphics, methods, preprocessCore (>= 1.29.0), RSQLite (>= 1.0.0), splines, stats, stats4, utils, zlibbioc

Enhances doMC, doMPI

LinkingTo preprocessCore

Suggests BSgenome.Hsapiens.UCSC.hg18, hapmap100kxba, pd.hg.u95av2, pd.mapping50k.xba240, pd.huex.1.0.st.v2, pd.hg18.60mer.expr, pd.hugene.1.0.st.v1, maqcExpression4plex, genefilter, limma, RColorBrewer, oligoData, BiocStyle, knitr, RUnit, biomaRt, AnnotationDbi, ACME, RCurl

VignetteBuilder knitr

Description A package to analyze oligonucleotide arrays (expression/SNP/tiling/exon) at probe-level. It currently supports Affymetrix (CEL files) and NimbleGen arrays (XYS files).

License LGPL (>= 2)

Collate AllGenerics.R methods-GenericArrays.R methods-GeneFeatureSet.R methods-ExpressionSet.R methods-ExpressionFeatureSet.R methods-ExpressionSet.R methods-LDS.R methods-FeatureSet.R methods-SnpFeatureSet.R methods-SnpCnvFeatureSet.R methods-TilingFeatureSet.R methods-HtaFeatureSet.R methods-DBPDInfo.R methods-background.R methods-normalization.R methods-summarization.R read.celfiles.R read.xysfiles.R utils-general.R utils-selectors.R todo-snp.R functions-crlmm.R

2 Contents

functions-snprma.R justSNPRMA.R justCRLMM.R methods-snp6.R
methods-genotype.R methods-PLMset.R zzz.R

LazyLoad Yes

biocViews Microarray, OneChannel, TwoChannel, Preprocessing, SNP, DifferentialExpression, ExonArray, GeneExpression, DataImport

git_url https://git.bioconductor.org/packages/oligo

git_branch RELEASE_3_19

git_last_commit f9b9ea0

git_last_commit_date 2024-04-30

Repository Bioconductor 3.19

Date/Publication 2024-05-21

Contents

oligo-package	3
basecontent	4
basicPLM	4
basicRMA	6
boxplot	7
chromosome	8
crlmm	8
darkColors	9
fitProbeLevelModel	10
getAffinitySplineCoefficients	11
getBaseProfile	12
getContainer	12
getCrlmmSummaries	13
getNetAffx	13
6. 6	14
6	15
getProbeInfo	15
getX	16
hist	17
image	18
justSNPRMA	19
list.xysfiles	19
MAplot	20
mm	22
	23
mmSequence	24
oligo-defunct	24
oligoPLM-class	25
paCalls	27
plotM-methods	29
pmAllele	29

oligo-package	2
UH9U-DACKA96	,

	pmFragmentLength																				30
	pmPosition																				30
	pmStrand																				31
	probeNames																				31
	read.celfiles																				32
	read.xysfiles																				33
	readSummaries																				35
	rma-methods																				35
	runDate																				37
	sequenceDesignMat	trix																			38
	snprma																				38
	summarize																				39
Index																					41
olig	o-package	The olig	o pa	ıckaş	ge:	a to	ool .	for	lov	w-l€	eve	l ai	naly	vsis	of	ol	igo	nu	cle	otio	de

Description

The **oligo** package provides tools to preprocess different oligonucleotide arrays types: expression, tiling, SNP and exon chips. The supported manufacturers are Affymetrix and NimbleGen.

It offers support to large datasets (when the **bigmemory** is loaded) and can execute preprocessing tasks in parallel (if, in addition to **bigmemory**, the **snow** package is also loaded).

Details

The package will read the raw intensity files (CEL for Affymetrix; XYS for NimbleGen) and allow the user to perform analyses starting at the feature-level.

Reading in the intensity files require the existence of data packages that contain the chip specific information (X/Y coordinates; feature types; sequence). These data packages packages are built using the **pdInfoBuilder** package.

For Affymetrix SNP arrays, users are asked to download the already built annotation packages from BioConductor. This is because these packages contain metadata that are not automatically created. The following annotation packages are available:

50 K Xba - pd.mapping 50 kxba. 240 50 K Hind - pd.mapping 50 khind. 240 250 K Sty - pd.mapping 250 k.sty 250 K Nsp - pd.mapping 250 k.nsp Genome Wide Snp 5 (SNP 5.0) - pd.genome wide snp. 5 Genome Wide Snp 6 (SNP 6.0) - pd.genome wide snp. 6

For users interested in genotype calls for SNP 5.0 and 6.0 arrays, we strongly recommend the use use the **crlmm** package, which implements a more efficient version of CRLMM.

Author(s)

Benilton Carvalho - <carvalho@bclab.org>

4 basicPLM

References

Carvalho, B.; Bengtsson, H.; Speed, T. P. & Irizarry, R. A. Exploration, Normalization, and Genotype Calls of High Density Oligonucleotide SNP Array Data. Biostatistics, 2006.

basecontent

Sequence Base Contents

Description

Function to compute the amounts of each nucleotide in a sequence.

Usage

basecontent(seq)

Arguments

seq

character vector of length n containg a valid sequence (A/T/C/G)

Value

matrix with n rows and 4 columns with the counts for each base.

Examples

```
sequences <- c("ATATATCCCCG", "TTTCCGAGC")
basecontent(sequences)</pre>
```

basicPLM

Simplified interface to PLM.

Description

Simplified interface to PLM.

Usage

```
basicPLM(pmMat, pnVec, normalize = TRUE, background = TRUE, transfo =
  log2, method = c('plm', 'plmr', 'plmrr', 'plmrc'), verbose = TRUE)
```

basicPLM 5

Arguments

pmMat Matrix of intensities to be processed.

pnVec Probeset names

normalize Logical flag: normalize?

background Logical flag: background adjustment?

transfo function: function to be used for data transformation prior to summarization.

method Name of the method to be used for normalization. 'plm' is the usual PLM model;

'plmr' is the (row and column) robust version of PLM; 'plmrr' is the row-robust

version of PLM; 'plmrc' is the column-robust version of PLM.

verbose Logical flag: verbose.

Value

A list with the following components:

Estimates A (length(pnVec) x ncol(pmMat)) matrix with probeset summaries.

StdErrors A (length(pnVec) x ncol(pmMat)) matrix with standard errors of 'Estimates'.

Residuals A (nrow(pmMat) x ncol(pmMat)) matrix of residuals.

Note

Currently, only RMA-bg-correction and quantile normalization are allowed.

Author(s)

Benilton Carvalho

See Also

rcModelPLMr, rcModelPLMrr, rcModelPLMrr, rcModelPLMrc, basicRMA

Examples

```
set.seed(1)
pms <- 2^matrix(rnorm(1000), nc=20)
colnames(pms) <- paste("sample", 1:20, sep="")
pns <- rep(letters[1:10], each=5)
res <- basicPLM(pms, pns, TRUE, TRUE)
res[['Estimates']][1:4, 1:3]
res[['StdErrors']][1:4, 1:3]
res[['Residuals']][1:20, 1:3]</pre>
```

6 basicRMA

basicRMA	Simplified interface to RMA.

Description

Simple interface to RMA.

Usage

```
basicRMA(pmMat, pnVec, normalize = TRUE, background = TRUE, bgversion = 2, destructive = FALSE, verbose
```

Arguments

pmMat Matrix of intensities to be processed.

pnVec Probeset names.

normalize Logical flag: normalize?

background Logical flag: background adjustment?

bgversion Version of background correction.

destructive Logical flag: use destructive methods?

verbose Logical flag: verbose.
... Not currently used.

Value

Matrix.

Examples

```
set.seed(1)
pms <- 2^matrix(rnorm(1000), nc=20)
colnames(pms) <- paste("sample", 1:20, sep="")
pns <- rep(letters[1:10], each=5)
res <- basicRMA(pms, pns, TRUE, TRUE)
res[, 1:3]</pre>
```

boxplot 7

Description

Boxplot for observed (log-)intensities in a FeatureSet-like object (ExpressionFeatureSet, ExonFeatureSet, SnpFeatureSet, TilingFeatureSet) and ExpressionSet.

Usage

```
## S4 method for signature 'FeatureSet'
boxplot(x, which=c("pm", "mm", "bg", "both",
"all"), transfo=log2, nsample=10000, target = "mps1", ...)
## S4 method for signature 'ExpressionSet'
boxplot(x, which, transfo=identity, nsample=10000, ...)
```

Arguments

X	a FeatureSet-like object or ExpressionSet object.
which	character defining what probe types are to be used in the plot.
transfo	a function to transform the data before plotting. See 'Details'
nsample	number of units to sample and build the plot.
	arguments to be passed to the default boxplot method.

Details

The 'transfo' argument will set the transformation to be used. For raw data, 'transfo=log2' is a common practice. For summarized data (which are often in log2-scale), no transformation is needed (therefore 'transfo=identity').

Note

The boxplot methods for FeatureSet and Expression use a sample (via sample) of the probes/probesets to produce the plot. Therefore, the user interested in reproducibility is advised to use set.seed.

See Also

```
hist, image, sample, set.seed
```

8 crlmm

chromosome

Accessor for chromosome information

Description

Returns chromosome information.

Usage

```
pmChr(object)
```

Arguments

object

TilingFeatureSet or SnpCallSet object

Details

chromosome() returns the chromosomal information for all probes and pmChr() subsets the output to the PM probes only (if a TilingFeatureSet object).

Value

Vector with chromosome information.

crlmm

Genotype Calls

Description

Performs genotype calls via CRLMM (Corrected Robust Linear Model with Maximum-likelihood based distances).

Usage

darkColors 9

Arguments

filenames character vector with the filenames.

outdir directory where the output (and some tmp files) files will be saved.
batch_size integer defining how many SNPs should be processed at a time.

recalibrate Logical - should recalibration be performed?

balance Control parameter to balance homozygotes and heterozygotes calls.

minLLRforCalls Minimum thresholds for genotype calls.

verbose Logical.

phenoData phenoData object or NULL
pkgname alt. pdInfo package to be used
reference logical, defaulting to TRUE ...

tmpdir Directory where temporary files are going to be stored at.

Value

SnpCallSetPlus object.

darkColors

Create set of colors, interpolating through a set of preferred colors.

Description

Create set of colors, interpolating through a set of preferred colors.

Usage

darkColors(n)
seqColors(n)
seqColors2(n)
divColors(n)

Arguments

n integer determining number of colors to be generated

Details

darkColors is based on the Dark2 palette in RColorBrewer, therefore useful to describe qualitative features of the data.

seqColors is based on Blues and generates a gradient of blues, therefore useful to describe quantitative features of the data. seqColors2 behaves similarly, but it is based on OrRd (white-orange-red). divColors is based on the RdBu pallete in RColorBrewer, therefore useful to describe quantitative

features ranging on two extremes.

10 fitProbeLevelModel

Examples

```
x <- 1:10
y <- 1:10
cols1 <- darkColors(10)
cols2 <- seqColors(10)
cols3 <- divColors(10)
cols4 <- seqColors2(10)
plot(x, y, col=cols1, xlim=c(1, 13), pch=19, cex=3)
points(x+1, y, col=cols2, pch=19, cex=3)
points(x+2, y, col=cols3, pch=19, cex=3)
points(x+3, y, col=cols4, pch=19, cex=3)
abline(0, 1, lty=2)
abline(-1, 1, lty=2)
abline(-2, 1, lty=2)
abline(-3, 1, lty=2)</pre>
```

fitProbeLevelModel

Tool to fit Probe Level Models.

Description

Fits robust Probe Level linear Models to all the (meta)probesets in an FeatureSet. This is carried out on a (meta)probeset by (meta)probeset basis.

Usage

fitProbeLevelModel(object, background=TRUE, normalize=TRUE, target="core", method="plm", verbose=TRUE

Arguments

object FeatureSet object. background Do background correction? normalize Do normalization? character vector describing the summarization target. Valid values are: 'probetarget set', 'core' (Gene/Exon), 'full' (Exon), 'extended' (Exon). method summarization method to be used. verbose verbosity flag. return final value as an S4 object (oligoPLM) if TRUE. If FALSE, final value is **S4** returned as a list. subset to be passed down to getProbeInfo for subsetting. See subset for details.

Value

fitProbeLevelModel returns an oligoPLM object, if S4=TRUE; otherwise, it will return a list.

Note

This is the initial port of fitPLM to oligo. Some features found on the original work by Ben Bolstad (in the affyPLM package) may not be yet available. If you found one of this missing characteristics, please contact Benilton Carvalho.

Author(s)

This is a simplified port from Ben Bolstad's work implemented in the affyPLM package. Problems with the implementation in oligo should be reported to Benilton Carvalho.

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California, Berkeley.

See Also

```
rma, summarizationMethods, subset
```

Examples

```
if (require(oligoData)){
  data(nimbleExpressionFS)
  fit <- fitProbeLevelModel(nimbleExpressionFS)
  image(fit)
  NUSE(fit)
  RLE(fit)
}</pre>
```

```
getAffinitySplineCoefficients
```

Estimate affinity coefficients.

Description

Estimate affinity coefficients using sequence information and splines.

Usage

```
getAffinitySplineCoefficients(intensities, sequences)
```

Arguments

```
intensities Intensity matrix sequences Probe sequences
```

12 getContainer

Value

Matrix with estimated coefficients.

See Also

getBaseProfile

getBaseProfile

Compute and plot nucleotide profile.

Description

Computes and, optionally, lots nucleotide profile, describing the sequence effect on intensities.

Usage

```
getBaseProfile(coefs, probeLength = 25, plot = FALSE, ...)
```

Arguments

coefs affinity spline coefficients.

probeLength length of probes

plot logical. Plots profile?

... arguments to be passed to matplot.

Value

Invisibly returns a matrix with estimated effects.

getContainer

Get container information for NimbleGen Tiling Arrays.

Description

Get container information for NimbleGen Tiling Arrays. This is useful for better identification of control probes.

Usage

```
getContainer(object, probeType)
```

Arguments

object A TilingFeatureSet or TilingFeatureSet object.

probeType String describing which probes to query ('pm', 'bg')

getCrlmmSummaries 13

Value

'character' vector with container information.

getCrlmmSummaries

Function to get CRLMM summaries saved to disk

Description

This will read the summaries written to disk and return them to the user as a SnpCallSetPlus or SnpCnvCallSetPlus object.

Usage

```
getCrlmmSummaries(tmpdir)
```

Arguments

tmpdir

directory where CRLMM saved the results to.

Value

If the data were from SNP 5.0 or 6.0 arrays, the function will return a SnpCnvCallSetPlus object. It will return a SnpCallSetPlus object, otherwise.

getNetAffx

NetAffx Biological Annotations

Description

Gets NetAffx Biological Annotations saved in the annotation package (Exon and Gene ST Affymetrix arrays).

Usage

```
getNetAffx(object, type = "probeset")
```

Arguments

object 'ExpressionSet' object (eg., result of rma())

type Either 'probeset' or 'transcript', depending on what type of summaries were

obtained.

14 getNgsColorsInfo

Details

This retrieves NetAffx annotation saved in the (pd) annotation package - annotation(object). It is only available for Exon ST and Gene ST arrays.

The 'type' argument should match the summarization target used to generate 'object'. The 'rma' method allows for two targets: 'probeset' (target='probeset') and 'transcript' (target='core', target='full', target='extended').

Value

'AnnotatedDataFrame' that can be used as featureData(object)

Author(s)

Benilton Carvalho

getNgsColorsInfo Helper function to extract color information for filenames on Nimble-Gen arrays.

Description

This function will (try to) extract the color information for NimbleGen arrays. This is useful when using read.xysfiles2 to parse XYS files for Tiling applications.

Usage

```
getNgsColorsInfo(path = ".", pattern1 = "_532", pattern2 = "_635", ...)
```

Arguments

path path where to look for files

pattern1 pattern to match files supposed to go to the first channel

pattern2 pattern to match files supposed to go to the second channel

extra arguments for list.xysfiles

Details

Many NimbleGen samples are identified following the pattern sampleID_532.XYS / sampleID_635.XYS. The function suggests sample names if all the filenames follow the standard above.

Value

A data.frame with, at least, two columns: 'channel1' and 'channel2'. A third column, 'sample-Names', is returned if the filenames follow the sampleID_532.XYS / sampleID_635.XYS standard.

Author(s)

Benilton Carvalho <bcarvalh@jhsph.edu>

getPlatformDesign 15

+D1-+	D -4-: Dl -4f D:1:4
getPlatformDesign	Retrieve Platform Design object

Description

Retrieve platform design object.

Usage

```
getPlatformDesign(object)
getPD(object)
```

Arguments

object FeatureSet object

Details

Retrieve platform design object.

Value

platformDesign or PDInfo object.

getProbeInfo	Probe information selector.

Description

A tool to simplify the selection of probe information, so user does not need to use the SQL approaches.

Usage

```
getProbeInfo(object, field, probeType = "pm", target = "core", sortBy = c("fid", "man_fsetid", "none"),
```

Arguments

object	FeatureSet object.
field	character string with names of $field(s)$ of interest to be obtained from database.
probeType	character string: 'pm' or 'mm'
target	Used only for Exon or Gene ST arrays: 'core', 'full', 'extended', 'probeset'.
sortBy	Field to be used for sorting.
	Arguments to be passed to subset

getX

Value

A data.frame with the probe level information.

Note

The code allows for querying info on MM probes, however it has been used mostly on PM probes.

Author(s)

Benilton Carvalho

Examples

```
if (require(oligoData)){
   data(affyGeneFS)
   availProbeInfo(affyGeneFS)
   probeInfo <- getProbeInfo(affyGeneFS, c('fid', 'x', 'y', 'chrom'))
   head(probeInfo)
   ## Selecting antigenomic background probes
   agenGene <- getProbeInfo(affyGeneFS, field=c('fid', 'fsetid', 'type'), target='probeset', subset= type == 'contr
   head(agenGene)
}</pre>
```

getX

Accessors for physical array coordinates.

Description

Accessors for physical array coordinates.

Usage

```
getX(object, type)
getY(object, type)
```

Arguments

object FeatureSet object

type 'character' defining the type of the probes to be queried. Valid options are 'pm',
 'mm', 'bg'

Value

A vector with the requested coordinates.

hist 17

Examples

```
## Not run:
x <- read.celfiles(list.celfiles())
theXpm <- getX(x, "pm")
theYpm <- getY(x, "pm")
## End(Not run)</pre>
```

hist

Density estimate

Description

Plot the density estimates for each sample

Usage

Arguments

Χ	FeatureSet or ExpressionSet object
transfo	a function to transform the data before plotting. See 'Details'.
nsample	number of units to sample and build the plot.
which	set of probes to be plotted ("pm", "mm", "bg", "both", "all").
	arguments to be passed to matplot

Details

The 'transfo' argument will set the transformation to be used. For raw data, 'transfo=log2' is a common practice. For summarized data (which are often in log2-scale), no transformation is needed (therefore 'transfo=identity').

Note

The hist methods for FeatureSet and Expression use a sample (via sample) of the probes/probesets to produce the plot (unless nsample > nrow(x)). Therefore, the user interested in reproducibility is advised to use set.seed.

18 image

image

Display a pseudo-image of a microarray chip

Description

Produces a pseudo-image (graphics::image) for each sample.

Usage

Arguments

FeatureSet object Х which integer indices of samples to be plotted (optional). transfo function to be applied to the data prior to plotting. Type of statistics to be used. type use.log Use log. add.legend Add legend. standardize Standardize residuals. Colors to be used. col Main title. main parameters to be passed to image

Examples

```
if(require(oligoData) & require(pd.hg18.60mer.expr)){
  data(nimbleExpressionFS)
  par(mfrow=c(1, 2))
  image(nimbleExpressionFS, which=4)
## fit <- fitPLM(nimbleExpressionFS)
## image(fit, which=4)
  plot(1) ## while fixing fitPLM TODO
}</pre>
```

justSNPRMA 19

justSNPRMA

Summarization of SNP data

Description

This function implements the SNPRMA method for summarization of SNP data. It works directly with the CEL files, saving memory.

Usage

```
justSNPRMA(filenames, verbose = TRUE, phenoData = NULL, normalizeToHapmap = TRUE)
```

Arguments

filenames character vector with the filenames.

verbose logical flag for verbosity.
phenoData a phenoData object or NULL

normalizeToHapmap

Normalize to Hapmap? Should always be TRUE, but it's kept here for future

use.

Value

SnpQSet or a SnpCnvQSet, depending on the array type.

Examples

```
## snprmaResults <- justSNPRMA(list.celfiles())</pre>
```

list.xysfiles

List XYS files

Description

Lists the XYS files.

Usage

```
list.xysfiles(...)
```

Arguments

... parameters to be passed to list.files

Details

The functions interface list.files and the user is asked to check that function for further details.

Value

Character vector with the filenames.

See Also

```
list.files
```

Examples

```
list.xysfiles()
```

MAplot

MA plots

Description

Create MA plots using a reference array (if one channel) or using channel2 as reference (if two channel).

Usage

```
MAplot(object, ...)
## S4 method for signature 'FeatureSet'
MAplot(object, what=pm, transfo=log2, groups,
       refSamples, which, pch=".", summaryFun=rowMedians,
       plotFun=smoothScatter, main="vs pseudo-median reference chip",
       pairs=FALSE, ...)
## S4 method for signature 'TilingFeatureSet'
MAplot(object, what=pm, transfo=log2, groups,
       refSamples, which, pch=".", summaryFun=rowMedians,
       plotFun=smoothScatter, main="vs pseudo-median reference chip",
       pairs=FALSE, ...)
## S4 method for signature 'PLMset'
MAplot(object, what=coefs, transfo=identity, groups,
       refSamples, which, pch=".", summaryFun=rowMedians,
       plotFun=smoothScatter, main="vs pseudo-median reference chip",
       pairs=FALSE, ...)
## S4 method for signature 'matrix'
MAplot(object, what=identity, transfo=identity,
```

MAplot 21

```
groups, refSamples, which, pch=".", summaryFun=rowMedians,
plotFun=smoothScatter, main="vs pseudo-median reference chip",
pairs=FALSE, ...)

## S4 method for signature 'ExpressionSet'
MAplot(object, what=exprs, transfo=identity,
    groups, refSamples, which, pch=".", summaryFun=rowMedians,
    plotFun=smoothScatter, main="vs pseudo-median reference chip",
    pairs=FALSE, ...)
```

Arguments

object FeatureSet, PLMset or ExpressionSet object. what function to be applied on object that will extract the statistics of interest, from which log-ratios and average log-intensities will be computed. transfo function to transform the data prior to plotting. groups factor describing groups of samples that will be combined prior to plotting. If missing, MvA plots are done per sample. integers (indexing samples) to define which subjects will be used to compute the refSamples reference set. If missing, a pseudo-reference chip is estimated using summaryFun. which integer (indexing samples) describing which samples are to be plotted. pch same as pch in plot summaryFun function that operates on a matrix and returns a vector that will be used to summarize data belonging to the same group (or reference) on the computation of grouped-stats. plotFun function to be used for plotting. Usually smoothScatter, plot or points. main string to be used in title. logical flag to determine if a matrix of MvA plots is to be generated pairs

Details

MAplot will take the following extra arguments:

1. subset: indices of elements to be plotted to reduce impact of plotting 100's thousands points (if pairs=FALSE only);

Other arguments to be passed downstream, like plot arguments.

- 2. span: see loess;
- 3. family.loess: see loess;
- 4. addLoess: logical flag (default TRUE) to add a loess estimate;
- 5. parParams: list of params to be passed to par() (if pairs=TRUE only);

Value

Plot

22 mm

Author(s)

Benilton Carvalho - based on Ben Bolstad's original MAplot function.

See Also

```
plot, smoothScatter
```

Examples

```
if(require(oligoData) & require(pd.hg18.60mer.expr)){
  data(nimbleExpressionFS)
  nimbleExpressionFS
  groups <- factor(rep(c('brain', 'UnivRef'), each=3))
  data.frame(sampleNames(nimbleExpressionFS), groups)
  MAplot(nimbleExpressionFS, pairs=TRUE, ylim=c(-.5, .5), groups=groups)
}</pre>
```

mm

Accessors and replacement methods for the intensity/PM/MM/BG matrices.

Description

Accessors and replacement methods for the PM/MM/BG matrices.

Usage

```
intensity(object)
mm(object, subset = NULL, target='core')
pm(object, subset = NULL, target='core')
bg(object, subset = NULL)
mm(object, subset = NULL, target='core')<-value
pm(object, subset = NULL, target='core')<-value
bg(object)<-value</pre>
```

Arguments

```
object FeatureSet object.

subset Not implemented yet.

value matrix object.

target One of 'probeset', 'core', 'full', 'extended'. This is ignored if the array design is something other than Gene ST or Exon ST.
```

mmindex 23

Details

For all objects but TilingFeatureSet, these methods will return matrices. In case of TilingFeatureSet objects, the value is a 3-dimensional array (probes x samples x channels).

intensity will return the whole intensity matrix associated to the object. pm, mm, bg will return the respective PM/MM/BG matrix.

When applied to ExonFeatureSet or GeneFeatureSet objects, pm will return the PM matrix at the transcript level ('core' probes) by default. The user should set the target argument accordingly if something else is desired. The valid values are: 'probeset' (Exon and Gene arrays), 'core' (Exon and Gene arrays), 'full' (Exon arrays) and 'extended' (Exon arrays).

The target argument has no effects when used on designs other than Gene and Exon ST.

Examples

```
if (require(maqcExpression4plex) & require(pd.hg18.60mer.expr)){
  xysPath <- system.file("extdata", package="maqcExpression4plex")
  xysFiles <- list.xysfiles(xysPath, full.name=TRUE)
  ngsExpressionFeatureSet <- read.xysfiles(xysFiles)
  pm(ngsExpressionFeatureSet)[1:10,]
}</pre>
```

mmindex

Accessors for PM, MM or background probes indices.

Description

Extracts the indexes for PM, MM or background probes.

Usage

```
mmindex(object, ...)
pmindex(object, ...)
bgindex(object, ...)
```

Arguments

```
object FeatureSet or DBPDInfo object
... Extra arguments, not yet implemented
```

Details

The indices are ordered by 'fid', i.e. they follow the order that the probes appear in the CEL/XYS files

Value

A vector of integers representing the rows of the intensity matrix that correspond to PM, MM or background probes.

24 oligo-defunct

Examples

```
## How pm() works
## Not run:
x <- read.celfiles(list.celfiles())
pms0 <- pm(x)
pmi <- pmindex(x)
pms1 <- exprs(x)[pmi,]
identical(pms0, pms1)
## End(Not run)</pre>
```

 ${\it mmSequence}$

Probe Sequeces

Description

Accessor to the (PM/MM/background) probe sequences.

Usage

```
mmSequence(object)
pmSequence(object, ...)
bgSequence(object, ...)
```

Arguments

object FeatureSet, AffySNPPDInfo or DBPDInfo object additional arguments

Value

A DNAStringSet containing the PM/MM/background probe sequence associated to the array.

oligo-defunct

Defunct Functions in Package 'oligo'

Description

The functions or variables listed here are no longer part of 'oligo'

Usage

```
fitPLM(...)
coefs(...)
resids(...)
```

oligoPLM-class 25

Arguments

... Arguments.

Details

fitPLM was replaced by fitProbeLevelModel, allowing faster execution and providing more specific models. fitPLM was based in the code written by Ben Bolstad in the affyPLM package. However, all the model-fitting functions are now in the package preprocessCore, on which fitProbeLevelModel depends.

coefs and resids, like fitPLM, were inherited from the affyPLM package. They were replaced respectively by coef and residuals, because this is how these statistics are called everywhere else in R.

oligoPLM-class

Class "oligoPLM"

Description

A class to represent Probe Level Models.

Objects from the Class

Objects can be created by calls of the form fitProbeLevelModel(FeatureSetObject), where FeatureSetObject is an object obtained through read.celfiles or read.xysfiles, representing intensities observed for different probes (which are grouped in probesets or meta-probesets) across distinct samples.

Slots

```
chip.coefs: "matrix" with chip/sample effects - probeset-level description: "MIAME" compliant description information.
phenoData: "AnnotatedDataFrame" with phenotypic data.
protocolData: "AnnotatedDataFrame" with protocol data.
probe.coefs: "numeric" vector with probe effects
weights: "matrix" with weights - probe-level
residuals: "matrix" with residuals - probe-level
se.chip.coefs: "matrix" with standard errors for chip/sample coefficients
se.probe.coefs: "numeric" vector with standard errors for probe effects
residualSE: scale - residual standard error
geometry: array geometry used for plots
method: "character" string describing method used for PLM
manufacturer: "character" string with manufacturer name
```

26 oligoPLM-class

```
annotation: "character" string with the name of the annotation package
narrays: "integer" describing the number of arrays
nprobes: "integer" describing the number of probes before summarization
nprobesets: "integer" describing the number of probesets after summarization
```

Methods

```
annotation signature(object = "oligoPLM"): accessor/replacement method to annotation slot
boxplot signature(x = "oligoPLM"): boxplot method
coef signature(object = "oligoPLM"): accessor/replacement method to coef slot
coefs.probe signature(object = "oligoPLM"): accessor/replacement method to coefs.probe slot
geometry signature(object = "oligoPLM"): accessor/replacement method to geometry slot
image signature(x = "oligoPLM"): image method
manufacturer signature(object = "oligoPLM"): accessor/replacement method to manufacturer
    slot
method signature(object = "oligoPLM"): accessor/replacement method to method slot
ncol signature(x = "oligoPLM"): accessor/replacement method to ncol slot
nprobes signature(object = "oligoPLM"): accessor/replacement method to nprobes slot
nprobesets signature(object = "oligoPLM"): accessor/replacement method to nprobesets slot
residuals signature(object = "oligoPLM"): accessor/replacement method to residuals slot
residualSE signature(object = "oligoPLM"): accessor/replacement method to residualSE slot
se signature(object = "oligoPLM"): accessor/replacement method to se slot
se.probe signature(object = "oligoPLM"): accessor/replacement method to se.probe slot
show signature(object = "oligoPLM"): show method
weights signature(object = "oligoPLM"): accessor/replacement method to weights slot
NUSE signature(x = "oligoPLM"): Boxplot of Normalized Unscaled Standard Errors (NUSE)
    or NUSE values.
RLE signature(x = "oligoPLM"): Relative Log Expression boxplot or values.
opset2eset signature(x = "oligoPLM") : Convert to ExpressionSet.
```

Author(s)

This is a port from Ben Bolstad's work implemented in the affyPLM package. Problems with the implementation in oligo should be reported to the package's maintainer.

References

Bolstad, BM (2004) Low Level Analysis of High-density Oligonucleotide Array Data: Background, Normalization and Summarization. PhD Dissertation. University of California, Berkeley.

See Also

rma, summarize

27 paCalls

Examples

```
## TODO: review code and fix broken
## Not run:
if (require(oligoData)){
  data(nimbleExpressionFS)
  fit <- fitProbeLevelModel(nimbleExpressionFS)</pre>
  image(fit)
  NUSE(fit)
  RLE(fit)
}
## End(Not run)
```

paCalls

Methods for P/A Calls

Description

Methods for Present/Absent Calls are meant to provide means of assessing whether or not each of the (PM) intensities are compatible with observations generated by background probes.

Usage

```
paCalls(object, method, ..., verbose=TRUE)
## S4 method for signature 'ExonFeatureSet'
paCalls(object, method, verbose = TRUE)
## S4 method for signature 'GeneFeatureSet'
paCalls(object, method, verbose = TRUE)
## S4 method for signature 'ExpressionFeatureSet'
paCalls(object, method, ..., verbose = TRUE)
```

Arguments

object Exon/Gene/Expression-FeatureSet object. String defining what method to use. See 'Details'. method Additional arguments passed to MAS5. See 'Details' verbose Logical flag for verbosity.

Details

For Whole Transcript arrays (Exon/Gene) the valid options for method are 'DABG' (p-values for each probe) and 'PSDABG' (p-values for each probeset). For Expression arrays, the only option currently available for method is 'MAS5'.

ABOUT MAS5 CALLS:

The additional arguments that can be passed to MAS5 are:

1. alpha1: a significance threshold in (0, alpha2);

28 paCalls

- 2. alpha2: a significance threshold in (alpha1, 0.5);
- 3. tau: a small positive constant;
- 4. ignore.saturated: if TRUE, do the saturation correction described in the paper, with a saturation level of 46000;

This function performs the hypothesis test:

H0: median(Ri) = tau, corresponding to absence of transcript H1: median(Ri) > tau, corresponding to presence of transcript

where Ri = (PMi - MMi) / (PMi + MMi) for each i a probe-pair in the probe-set represented by data.

The p-value that is returned estimates the usual quantity:

Pr(observing a more "present looking" probe-set than data | data is absent)

So that small p-values imply presence while large ones imply absence of transcript. The detection call is computed by thresholding the p-value as in:

call "P" if p-value < alpha1 call "M" if alpha1 <= p-value < alpha2 call "A" if alpha2 <= p-value

Value

A matrix (of dimension dim(PM) if method="DABG" or "MAS5"; of dimension length(unique(probeNames(object))) x ncol(object) if method="PSDABG") with p-values for P/A Calls.

Author(s)

Benilton Carvalho

References

Clark et al. Discovery of tissue-specific exons using comprehensive human exon microarrays. Genome Biol (2007) vol. 8 (4) pp. R64

Liu, W. M. and Mei, R. and Di, X. and Ryder, T. B. and Hubbell, E. and Dee, S. and Webster, T. A. and Harrington, C. A. and Ho, M. H. and Baid, J. and Smeekens, S. P. (2002) Analysis of high density expression microarrays with signed-rank call algorithms, Bioinformatics, 18(12), pp. 1593–1599.

Liu, W. and Mei, R. and Bartell, D. M. and Di, X. and Webster, T. A. and Ryder, T. (2001) Rank-based algorithms for analysis of microarrays, Proceedings of SPIE, Microarrays: Optical Technologies and Informatics, 4266.

Affymetrix (2002) Statistical Algorithms Description Document, Affymetrix Inc., Santa Clara, CA, whitepaper. http://www.affymetrix.com/support/technical/whitepapers/sadd_whitepaper.pdf

Examples

```
## Not run:
if (require(oligoData) & require(pd.huex.1.0.st.v2)){
  data(affyExonFS)
  ## Get only 2 samples for example
  dabgP = paCalls(affyExonFS[, 1:2])
  dabgPS = paCalls(affyExonFS[, 1:2], "PSDABG")
```

plotM-methods 29

```
head(dabgP) ## for probe
head(dabgPS) ## for probeset
}
## End(Not run)
```

plotM-methods

Methods for Log-Ratio plotting

Description

The plotM methods are meant to plot log-ratios for different classes of data.

Methods

```
object = "SnpQSet", i = "character" Plot log-ratio for SNP data for sample i.
object = "SnpQSet", i = "integer" Plot log-ratio for SNP data for sample i.
object = "SnpQSet", i = "numeric" Plot log-ratio for SNP data for sample i.
object = "TilingQSet", i = "missing" Plot log-ratio for Tiling data for sample i.
```

pmAllele

Access the allele information for PM probes.

Description

Accessor to the allelic information for PM probes.

Usage

```
pmAllele(object)
```

Arguments

object

SnpFeatureSet or PDInfo object.

30 pmPosition

pmFragmentLength	Access the fragment length for PM probes.

Description

Accessor to the fragment length for PM probes.

Usage

```
pmFragmentLength(object, enzyme, type=c('snp', 'cn'))
```

Arguments

object PDInfo or SnpFeatureSet object.

enzyme Enzyme to be used for query. If missing, all enzymes are used.

type Type of probes to be used: 'snp' for SNP probes; 'cn' for Copy Number probes.

Value

A list of length equal to the number of enzymes used for digestion. Each element of the list is a data.frame containing:

- row: the row used to link to the PM matrix;
- length: expected fragment length.

Note

There is not a 1:1 relationship between probes and expected fragment length. For one enzyme, a given probe may be associated to multiple fragment lengths. Therefore, the number of rows in the data.frame may not match the number of PM probes and the row column should be used to match the fragment length with the PM matrix.

pmPosition Accessor to position information

Description

pmPosition will return the genomic position for the (PM) probes.

Usage

```
pmPosition(object)
pmOffset(object)
```

pmStrand 31

Arguments

object

AffySNPPDInfo, TilingFeatureSet or SnpCallSet object

Details

pmPosition will return genomic position for PM probes on a tiling array. pmOffset will return the offset information for PM probes on SNP arrays.

pmStrand

Accessor to the strand information

Description

Returns the strand information for PM probes (0 - sense / 1 - antisense).

Usage

```
pmStrand(object)
```

Arguments

object

AffySNPPDInfo or TilingFeatureSet object

probeNames

Accessor to feature names

Description

Accessors to featureset names.

Usage

```
probeNames(object, subset = NULL, ...)
probesetNames(object, ...)
```

Arguments

FeatureSet or DBPDInfo object not implemented yet. subset

Arguments (like 'target') passed to downstream methods. . . .

Value

probeNames returns a string with the probeset names for *each probe* on the array. probesetNames, on the other hand, returns the *unique probeset names*.

32 read.celfiles

read.celfiles	Parser to CEL files
---------------	---------------------

Description

Reads CEL files.

Usage

```
read.celfiles(..., filenames, pkgname, phenoData, featureData, experimentData, protocolData, notes, verbose=TRUE, sampleNames, rm.mask=FALSE, rm.outliers=FALSE, rm.extra=FALSE, checkType=TRUE) read.celfiles2(channel1, channel2, pkgname, phenoData, featureData, experimentData, protocolData, notes, verbose=TRUE, sampleNames, rm.mask=FALSE, rm.outliers=FALSE, rm.extra=FALSE, checkType=TRUE)
```

Arguments

... names of files to be read.

filenames a character vector with the CEL filenames.

channel1 a character vector with the CEL filenames for the first 'channel' on a Tiling

application

channel2 a character vector with the CEL filenames for the second 'channel' on a Tiling

application

pkgname alternative data package to be loaded.

phenoData phenoData
featureData featureData
experimentData experimentData
protocolData protocolData

notes notes verbose logical

sampleNames character vector with sample names (usually better descriptors than the file-

names)

rm.mask logical. Read masked?
rm.outliers logical. Remove outliers?
rm.extra logical. Remove extra?

checkType logical. Check type of each file? This can be time consuming.

read.xysfiles 33

Details

When using 'affyio' to read in CEL files, the user can read compressed CEL files (CEL.gz). Additionally, 'affyio' is much faster than 'affxparser'.

The function guesses which annotation package to use from the header of the CEL file. The user can also provide the name of the annotation package to be used (via the pkgname argument). If the annotation package cannot be loaded, the function returns an error. If the annotation package is not available from BioConductor, one can use the pdInfoBuilder package to build one.

Value

```
ExpressionFeatureSet
if Expresssion arrays
ExonFeatureSet if Exon arrays
SnpFeatureSet if SNP arrays
TilingFeatureSet
if Tiling arrays
```

See Also

```
list.celfiles, read.xysfiles
```

Examples

```
if(require(pd.mapping50k.xba240) & require(hapmap100kxba)){
celPath <- system.file("celFiles", package="hapmap100kxba")
celFiles <- list.celfiles(celPath, full.name=TRUE)
affySnpFeatureSet <- read.celfiles(celFiles)
}</pre>
```

read.xysfiles

Parser to XYS files

Description

NimbleGen provides XYS files which are read by this function.

Usage

```
read.xysfiles(..., filenames, pkgname, phenoData, featureData,
experimentData, protocolData, notes, verbose=TRUE, sampleNames,
checkType=TRUE)

read.xysfiles2(channel1, channel2, pkgname, phenoData, featureData,
experimentData, protocolData, notes, verbose=TRUE, sampleNames,
checkType=TRUE)
```

34 read.xysfiles

Arguments

... file names

filenames character vector with filenames.

channel1 a character vector with the XYS filenames for the first 'channel' on a Tiling

application

channel a character vector with the XYS filenames for the second 'channel' on a Tiling

application

pkgname character vector with alternative PD Info package name

phenoData phenoData
featureData featureData
experimentData experimentData
protocolData protocolData

notes notes verbose verbose

sampleNames character vector with sample names (usually better descriptors than the file-

names)

checkType logical. Check type of each file? This can be time consuming.

Details

The function will read the XYS files provided by NimbleGen Systems and return an object of class FeatureSet.

The function guesses which annotation package to use from the header of the XYS file. The user can also provide the name of the annotation package to be used (via the pkgname argument). If the annotation package cannot be loaded, the function returns an error. If the annotation package is not available from BioConductor, one can use the pdInfoBuilder package to build one.

Value

See Also

```
list.xysfiles, read.celfiles
```

Examples

```
if (require(maqcExpression4plex) & require(pd.hg18.60mer.expr)){
xysPath <- system.file("extdata", package="maqcExpression4plex")
xysFiles <- list.xysfiles(xysPath, full.name=TRUE)
ngsExpressionFeatureSet <- read.xysfiles(xysFiles)
}</pre>
```

readSummaries 35

l by crlmm	Read summaries	readSummaries
------------	----------------	---------------

Description

This function read the different summaries generated by crlmm.

Usage

```
readSummaries(type, tmpdir)
```

Arguments

type of summary of character class: 'alleleA', 'alleleB', 'alleleA-sense', 'alleleA-

antisense', 'alleleB-sense', 'alleleB-antisense', 'calls', 'llr', 'conf'.

tmpdir directory containing the output saved by crlmm

Details

On the 50K and 250K arrays, given a SNP, there are probes on both strands (sense and antisense). For this reason, the options 'alleleA-sense', 'alleleA-antisense', 'alleleB-sense' and 'alleleB-antisense' should be used **only** with such arrays (XBA, HIND, NSP or STY).

On the SNP 5.0 and SNP 6.0 platforms, this distinction does not exist in terms of algorithm (note that the actual strand could be queried from the annotation package). For these arrays, options 'alleleA', 'alleleB' are the ones to be used.

The options calls, 11r and conf will return, respectively, the CRLMM calls, log-likelihood ratios (for devel purpose **only**) and CRLMM confidence calls matrices.

Value

Matrix with values of summaries.

Description

Robust Multichip Average preprocessing methodology. This strategy allows background subtraction, quantile normalization and summarization (via median-polish).

36 rma-methods

Usage

```
## S4 method for signature 'ExonFeatureSet'
rma(object, background=TRUE, normalize=TRUE, subset=NULL, target="core")
## S4 method for signature 'HTAFeatureSet'
rma(object, background=TRUE, normalize=TRUE, subset=NULL, target="core")
## S4 method for signature 'ExpressionFeatureSet'
rma(object, background=TRUE, normalize=TRUE, subset=NULL)
## S4 method for signature 'GeneFeatureSet'
rma(object, background=TRUE, normalize=TRUE, subset=NULL, target="core")
## S4 method for signature 'SnpCnvFeatureSet'
rma(object, background=TRUE, normalize=TRUE, subset=NULL)
```

Arguments

object Exon/HTA/Expression/Gene/SnpCnv-FeatureSet object.

background Logical - perform RMA background correction?

normalize Logical - perform quantile normalization?

subset To be implemented.

target Level of summarization (only for Exon/Gene arrays)

Methods

- signature(object = "ExonFeatureSet") When applied to an ExonFeatureSet object, rma can produce summaries at different levels: probeset (as defined in the PGF), core genes (as defined in the core.mps file), full genes (as defined in the full.mps file) or extended genes (as defined in the extended.mps file). To determine the level for summarization, use the target argument.
- signature(object = "ExpressionFeatureSet") When used on an ExpressionFeatureSet object, rma produces summaries at the probeset level (as defined in the CDF or NDF files, depending on the manufacturer).
- signature(object = "GeneFeatureSet") When applied to a GeneFeatureSet object, rma can produce summaries at different levels: probeset (as defined in the PGF) and 'core genes' (as defined in the core.mps file). To determine the level for summarization, use the target argument.
- signature(object = "HTAFeatureSet") When applied to a HTAFeatureSet object, rma can produce summaries at different levels: probeset (as defined in the PGF) and 'core genes' (as defined in the core.mps file). To determine the level for summarization, use the target argument.
- signature(object = "SnpCnvFeatureSet") If used on a SnpCnvFeatureSet object (ie., SNP 5.0 or SNP 6.0 arrays), rma will produce summaries for the CNV probes. Note that this is an experimental feature for internal (and quick) assessment of CNV probes. We recommend the use of the 'crlmm' package, which contains a Copy Number tool specifically designed for these data.

runDate 37

References

Rafael. A. Irizarry, Benjamin M. Bolstad, Francois Collin, Leslie M. Cope, Bridget Hobbs and Terence P. Speed (2003), Summaries of Affymetrix GeneChip probe level data Nucleic Acids Research 31(4):e15

Bolstad, B.M., Irizarry R. A., Astrand M., and Speed, T.P. (2003), A Comparison of Normalization Methods for High Density O ligonucleotide Array Data Based on Bias and Variance. Bioinformatics 19(2):185-193

Irizarry, RA, Hobbs, B, Collin, F, Beazer-Barclay, YD, Antonellis, KJ, Scherf, U, Speed, TP (2003) Exploration, Normalizati on, and Summaries of High Density Oligonucleotide Array Probe Level Data. Biostatistics. Vol. 4, Number 2: 249-264

See Also

snprma

Examples

```
if (require(maqcExpression4plex) & require(pd.hg18.60mer.expr)){
  xysPath <- system.file("extdata", package="maqcExpression4plex")
  xysFiles <- list.xysfiles(xysPath, full.name=TRUE)
  ngsExpressionFeatureSet <- read.xysfiles(xysFiles)
  summarized <- rma(ngsExpressionFeatureSet)
  show(summarized)
}</pre>
```

runDate

Date of scan

Description

Retrieves date information in CEL/XYS files.

Usage

```
runDate(object)
```

Arguments

object

'FeatureSet' object.

38 snprma

```
sequenceDesignMatrix Create design matrix for sequences
```

Description

Creates design matrix for sequences.

Usage

```
sequenceDesignMatrix(seqs)
```

Arguments

seqs

character vector of 25-mers.

Details

This assumes all sequences are 25bp long.

The design matrix is often used when the objective is to adjust intensities by sequence.

Value

Matrix with length(seqs) rows and 75 columns.

Examples

```
genSequence <- function(x)
    paste(sample(c("A", "T", "C", "G"), 25, rep=TRUE), collapse="", sep="")
seqs <- sapply(1:10, genSequence)
X <- sequenceDesignMatrix(seqs)
Y <- rnorm(10, mean=12, sd=2)
Ydemean <- Y-mean(Y)
X[1:10, 1:3]
fit <- lm(Ydemean~X)
coef(fit)</pre>
```

snprma

Preprocessing SNP Arrays

Description

This function preprocess SNP arrays.

Usage

```
snprma(object, verbose = TRUE, normalizeToHapmap = TRUE)
```

summarize 39

Arguments

object SnpFeatureSet object verbose Verbosity flag. logical

normalizeToHapmap internal

Value

A SnpQSet object.

summarize

Tools for microarray preprocessing.

Description

These are tools to preprocess microarray data. They include background correction, normalization and summarization methods.

Usage

```
backgroundCorrectionMethods()
normalizationMethods()
summarizationMethods()
backgroundCorrect(object, method=backgroundCorrectionMethods(), copy=TRUE, extra, subset=NULL, target
summarize(object, probes=rownames(object), method="medianpolish", verbose=TRUE, ...)
## S4 method for signature 'FeatureSet'
normalize(object, method=normalizationMethods(), copy=TRUE, subset=NULL, target='core', verbose=TRUE,
## S4 method for signature 'matrix'
normalize(object, method=normalizationMethods(), copy=TRUE, verbose=TRUE, ...)
## S4 method for signature 'ff_matrix'
normalize(object, method=normalizationMethods(), copy=TRUE, verbose=TRUE, ...)
normalizeToTarget(object, targetDist, method="quantile", copy=TRUE, verbose=TRUE)
```

Arguments

object Object containing probe intensities to be preprocessed.

method String determining which method to use at that preprocessing step.

targetDist Vector with the target distribution

probes Character vector that identifies the name of the probes represented by the rows

of object.

copy Logical flag determining if data must be copied before processing (TRUE), or if

data can be overwritten (FALSE).

subset Not yet implemented.

target One of the following values: 'core', 'full', 'extended', 'probeset'. Used only

with Gene ST and Exon ST designs.

40 summarize

extra Extra arguments to be passed to other methods.

verbose Logical flag for verbosity.

Arguments to be passed to methods.

Details

Number of rows of object must match the length of probes.

Value

backgroundCorrectionMethods and normalizationMethods will return a character vector with the methods implemented currently.

backgroundCorrect, normalize and normalizeToTarget will return a matrix with same dimensions as the input matrix. If they are applied to a FeatureSet object, the PM matrix will be used as input.

The summarize method will return a matrix with length(unique(probes)) rows and ncol(object) columns.

Examples

```
ns <- 100
nps <- 1000
np <- 10
intensities <- matrix(rnorm(ns*nps*np, 8000, 400), nc=ns)</pre>
ids <- rep(as.character(1:nps), each=np)</pre>
bgCorrected <- backgroundCorrect(intensities)</pre>
normalized <- normalize(bgCorrected)</pre>
summarizationMethods()
expression <- summarize(normalized, probes=ids)</pre>
intensities[1:20, 1:3]
expression[1:20, 1:3]
target <- rnorm(np*nps)</pre>
normalizedToTarget <- normalizeToTarget(intensities, target)</pre>
if (require(oligoData) & require(pd.hg18.60mer.expr)){
  ## Example of normalization with real data
  data(nimbleExpressionFS)
  boxplot(nimbleExpressionFS, main='Original')
  for (mtd in normalizationMethods()){
    message('Normalizing with ', mtd)
    res <- normalize(nimbleExpressionFS, method=mtd, verbose=FALSE)</pre>
    boxplot(res, main=mtd)
  }
}
```

Index

* IO	pmAllele, 29
read.celfiles, 32	pmFragmentLength, 30
read.xysfiles,33	pmPosition, 30
* classes	pmStrand, 31
oligoPLM-class, 25	probeNames, 31
* classif	readSummaries, 35
crlmm, 8	sequenceDesignMatrix, 38
getNetAffx, 13	snprma, 38
runDate, 37	summarize, 39
* file	* methods
list.xysfiles, 19	boxplot, 7
* hplot	hist, 17
boxplot, 7	MAplot, 20
darkColors,9	plotM-methods, 29
hist, 17	rma-methods, 35
image, 18	* package
MAplot, 20	oligo-package, 3
* loess	* smooth
MAplot, 20	MAplot, 20
* manip	
basecontent, 4	annotation, oligoPLM-method
basicPLM, 4	(oligoPLM-class), 25
basicRMA, 6	availProbeInfo (getProbeInfo), 15
chromosome, 8	backgroundCorrect(summarize), 39
fitProbeLevelModel, 10	backgroundCorrect,FeatureSet-method
<pre>getAffinitySplineCoefficients, 11</pre>	(summarize), 39
getBaseProfile, 12	backgroundCorrect,ff_matrix-method
getContainer, 12	(summarize), 39
getCrlmmSummaries, 13	backgroundCorrect, matrix-method
getNgsColorsInfo, 14	(summarize), 39
getPlatformDesign, 15	backgroundCorrect-methods (summarize),
getProbeInfo, 15	39
getX, 16	backgroundCorrectionMethods
justSNPRMA, 19	(summarize), 39
mm, 22	basecontent, 4
mmindex, 23	basicPLM, 4
mmSequence, 24	basicRMA, 5, 6
oligo-defunct, 24	bg (mm), 22
paCalls, 27	bg, FeatureSet-method (mm), 22
•	

INDEX

bg,TilingFeatureSet-method(mm),22	<pre>getAffinitySplineCoefficients, 11</pre>
bg<- (mm), 22	getBaseProfile, 12
<pre>bg<-,FeatureSet,ff_matrix-method(mm),</pre>	getContainer, 12
22	<pre>getContainer,TilingFeatureSet-method</pre>
bg<-,FeatureSet,matrix-method(mm),22	(getContainer), 12
bg<-,TilingFeatureSet,array-method	<pre>getContainer-methods (getContainer), 12</pre>
(mm), 22	getCrlmmSummaries, 13
bgindex (mmindex), 23	getNetAffx, 13
bgindex,DBPDInfo-method(mmindex),23	<pre>getNetAffx,ExpressionSet-method</pre>
bgindex, FeatureSet-method(mmindex), 23	(getNetAffx), 13
bgSequence (mmSequence), 24	<pre>getNetAffx-methods (getNetAffx), 13</pre>
bgSequence,DBPDInfo-method	<pre>getNgsColorsInfo, 14</pre>
(mmSequence), 24	<pre>getPD (getPlatformDesign), 15</pre>
bgSequence,ExonFeatureSet-method	getPlatformDesign, 15
(mmSequence), 24	<pre>getPlatformDesign,FeatureSet-method</pre>
bgSequence,FeatureSet-method	(getPlatformDesign), 15
(mmSequence), 24	getProbeInfo, 10, 15
bgSequence,GeneFeatureSet-method	getX, 16
(mmSequence), 24	<pre>getX,DBPDInfo-method(getX), 16</pre>
boxplot, 7	getX, FeatureSet-method(getX), 16
boxplot,ExpressionSet-method(boxplot),	getX-methods (getX), 16
7	getY (getX), 16
boxplot, FeatureSet-method (boxplot), 7	<pre>getY,DBPDInfo-method(getX), 16</pre>
boxplot,oligoPLM-method	getY, FeatureSet-method(getX), 16
(oligoPLM-class), 25	getY-methods (getX), 16
boxplot, PLMset-method (boxplot), 7	
boxplot-methods (boxplot), 7	hist, 7, 17
ch romocomo 0	hist,ExpressionSet-method(hist), 17
chromosome, 8	hist, FeatureSet-method (hist), 17
chromosome<- (chromosome), 8	hist-methods (hist), 17
chromosome<-, AnnotatedDataFrame, character-m	
(chromosome), 8	image, 7, 18
cleanPlatformName (read.celfiles), 32	image, FeatureSet-method (image), 18
<pre>coef, oligoPLM-method (oligoPLM-class),</pre>	<pre>image,oligoPLM-method(oligoPLM-class) 25</pre>
coefs(oligo-defunct), 24	image, PLMset-method (image), 18
coefs.probe(oligoPLM-class), 25	image-methods (image), 18
coefs.probe,oligoPLM-method	intensity (mm), 22
(oligoPLM-class), 25	intensity, FeatureSet-method (mm), 22
crlmm, 8	intensity<- (mm), 22
darkColors,9	<pre>intensity<-,FeatureSet-method (mm), 22</pre>
· · · · · · · · · · · · · · · · · · ·	• • • • • • • • • • • • • • • • • • • •
divColors(darkColors), 9	justCRLMM(crlmm), 8
FeatureSet, 10	justSNPRMA, 19
fitPLM(oligo-defunct), 24	
fitProbeLevelModel, 10	list.celfiles, 33
,	list.files, <i>19</i> , <i>20</i>
geometry,oligoPLM-method	list.xysfiles, 19, 34
(oligoPLM-class), 25	loess, <i>21</i>

INDEX 43

manufacturer,oligoPLM-method	nprobes (oligoPLM-class), 25
(oligoPLM-class), 25	nprobes,oligoPLM-method
MAplot, 20	(oligoPLM-class), 25
MAplot, ExpressionSet-method (MAplot), 20	nprobesets(oligoPLM-class), 25
MAplot, FeatureSet-method (MAplot), 20	nprobesets,oligoPLM-method
MAplot, matrix-method (MAplot), 20	(oligoPLM-class), 25
MAplot, PLMset-method (MAplot), 20	NUSE (oligoPLM-class), 25
MAplot,TilingFeatureSet-method	NUSE, oligoPLM-method (oligoPLM-class),
(MAplot), 20	25
MAplot-methods (MAplot), 20	
method(oligoPLM-class), 25	oligo-defunct, 24
method,oligoPLM-method	oligo-package, 3
(oligoPLM-class), 25	oligoPLM, 10
mm, 22	oligoPLM (oligoPLM-class), 25
mm, FeatureSet-method (mm), 22	oligoPLM-class, 25
mm, TilingFeatureSet-method (mm), 22	_
mm<- (mm), 22	opset2eset (oligoPLM-class), 25
mm<-,FeatureSet,ANY,ANY,ff_matrix-method	opset2eset,oligoPLM-method
(mm), 22	(oligoPLM-class), 25
mm<-,FeatureSet,ANY,ANY,matrix-method	
(mm), 22	paCalls, 27
mm<-,TilingFeatureSet,ANY,ANY,array-method	paCalls,ExonFeatureSet-method
(mm), 22	(paCalls), 27
mmindex, 23	paCalls,ExpressionFeatureSet-method
mmindex, DBPDInfo-method (mmindex), 23	(paCalls), 27
mmindex, FeatureSet-method (mmindex), 23	paCalls,GeneFeatureSet-method
mmSequence, 24	(paCalls), 27
mmSequence, AffySNPPDInfo-method	plot, 22
(mmSequence), 24	plotM(plotM-methods), 29
mmSequence, DBPDInfo-method	plotM,SnpQSet,character-method
(mmSequence), 24	(plotM-methods), 29
mmSequence, FeatureSet-method	plotM,SnpQSet,integer-method
(mmSequence), 24	(plotM-methods), 29
(plotM,SnpQSet,numeric-method
<pre>ncol,oligoPLM-method(oligoPLM-class),</pre>	(plotM-methods), 29
25	plotM,TilingQSet,missing-method
normalizationMethods (summarize), 39	(plotM-methods), 29
normalize, FeatureSet-method	plotM-methods, 29
(summarize), 39	pm (mm), 22
normalize, ff_matrix-method (summarize),	pm, FeatureSet-method (mm), 22
39	pm, GenericFeatureSet-method (mm), 22
normalize, matrix-method (summarize), 39	pm, TilingFeatureSet-method (mm), 22
normalizeToTarget (summarize), 39	pm<- (mm), 22
normalizeToTarget,ff_matrix-method	pm<-,FeatureSet,ANY,ANY,ff_matrix-method
(summarize), 39	(mm), 22
normalizeToTarget,matrix-method	pm<-,FeatureSet,ANY,ANY,matrix-method
(summarize), 39	(mm), 22
normalizeToTarget-methods (summarize),	pm<-,GenericFeatureSet,ANY,ANY,ff_matrix-method
39	(mm), 22

INDEX INDEX

pm<-, GenericFeatureSet, ANY, ANY, matrix-method	
(mm), 22	(mmSequence), 24
<pre>pm<-,TilingFeatureSet,ANY,ANY,array-method</pre>	<pre>pmSequence,FeatureSet-method</pre>
pmAllele, 29	pmSequence, GeneFeatureSet-method
pmAllele,AffySNPPDInfo-method	(mmSequence), 24
(pmAllele), 29	pmSequence, stArrayDBPDInfo-method
pmAllele, SnpFeatureSet-method	(mmSequence), 24
(pmAllele), 29	pmStrand, 31
	pmStrand, AffySNPPDInfo-method
pmChr (chromosome), 8	(pmStrand), 31
pmChr, ExonFeatureSet-method	pmStrand, TilingFeatureSet-method
(chromosome), 8	(pmStrand), 31
pmChr, FeatureSet-method (chromosome), 8	probeNames, 31
pmChr, GeneFeatureSet-method	
(chromosome), 8	probeNames, DBPDInfo-method
pmFragmentLength, 30	(probeNames), 31
pmFragmentLength,AffySNPPDInfo-method	probeNames, ExonFeatureSet-method
(pmFragmentLength), 30	(probeNames), 31
pmFragmentLength,SnpFeatureSet-method	probeNames, FeatureSet-method
(pmFragmentLength), 30	(probeNames), 31
pmindex (mmindex), 23	probeNames, GeneFeatureSet-method
<pre>pmindex, DBPDInfo-method (mmindex), 23</pre>	(probeNames), 31
pmindex, FeatureSet-method (mmindex), 23	probeNames, stArrayDBPDInfo-method
pmindex, GenericFeatureSet-method	(probeNames), 31
(mmindex), 23	probesetNames (probeNames), 31
<pre>pmindex, GenericPDInfo-method (mmindex),</pre>	probesetNames, FeatureSet-method
23	(probeNames), 31
pmindex,stArrayDBPDInfo-method	noModelDIM 5
(mmindex), 23	rcModelPLM, 5
pmOffset (pmPosition), 30	rcModelPLMr, 5
pmOffset,AffySNPPDInfo-method	rcModelPLMrc, 5
(pmPosition), 30	rcModelPLMrr, 5
pmPosition, 30	read.celfiles, 32, 34
pmPosition,ExpressionPDInfo-method	read.celfiles2(read.celfiles), 32
(pmPosition), 30	read.xysfiles, 33, 33
pmPosition, FeatureSet-method	read.xysfiles2(read.xysfiles), 33
	readSummaries, 35
(pmPosition), 30	resids (oligo-defunct), 24
pmPosition, TilingFeatureSet-method	residuals,oligoPLM-method
(pmPosition), 30	(oligoPLM-class), 25
pmPosition, TilingPDInfo-method	residualSE (oligoPLM-class), 25
(pmPosition), 30	residualSE,oligoPLM-method
pmSequence (mmSequence), 24	(oligoPLM-class), 25
pmSequence, AffyGenePDInfo-method	RLE (oligoPLM-class), 25
(mmSequence), 24	RLE, oligoPLM-method (oligoPLM-class), 25
pmSequence, AffySNPPDInfo-method	rma, 11, 26
(mmSequence), 24	rma (rma-methods), 35
pmSequence,DBPDInfo-method	rma,ExonFeatureSet-method
(mmSequence), 24	(rma-methods), 35

INDEX 45

```
rma, ExpressionFeatureSet-method
        (rma-methods), 35
rma, GeneFeatureSet-method
        (rma-methods), 35
rma, GenericFeatureSet-method
        (rma-methods), 35
rma, HTAFeatureSet-method (rma-methods),
        35
rma, SnpCnvFeatureSet-method
        (rma-methods), 35
rma-methods, 35
runDate, 37
runDate, FeatureSet-method (runDate), 37
runDate-methods (runDate), 37
sample, 7
se (oligoPLM-class), 25
se,oligoPLM-method(oligoPLM-class), 25
se.probe(oligoPLM-class), 25
se.probe,oligoPLM-method
        (oligoPLM-class), 25
seqColors (darkColors), 9
seqColors2 (darkColors), 9
sequenceDesignMatrix, 38
set.seed, 7
show,oligoPLM-method(oligoPLM-class),
smoothScatter, 22
snprma, 37, 38
subset, 10, 15
summarizationMethods, 11
\verb|summarizationMethods| (\verb|summarize|), 39|
summarize, 26, 39
summarize,ff_matrix-method(summarize),
summarize, matrix-method (summarize), 39
summarize-methods (summarize), 39
{\tt weights,oligoPLM-method}
        (oligoPLM-class), 25
```